



**Radiological Physics Associates**  
Medical Physics Consulting

2603 Elm Hill Pike, Suite P  
Nashville, TN 37214

800-262-6723  
inspections@rpanashville.com  
www.rpanashville.com

## Radiation Shielding Plan

January 2, 2026

**Kelly Judd**

**Medical Imaging Department**

Pikeville Urgent Care

3062 Main Street

Pikeville, Tennessee 37367

SUBJECT: Radiation Shielding Design for:

**Equipment:**

Toshiba Aquilion Prime

**Room:**

CT

**Dear Kelly,**

On the above indicated date, a radiation shielding plan was conducted for your facility according to the following standards:

- National Council on Radiation Protection and Measurements, NCRP Report No. 147 – Structural Shielding Design for Medical X-Ray Imaging Facilities

- American Association of Physicists in Medicine (AAPM) established protocol

Modifications to the specifications and calculations recommended by the organizations mentioned above are incorporated to comply with current Federal and State regulations.

The shielding design goals are such that the permissible exposure in controlled areas is no greater than 10 mR/week and no greater than 2 mR/week in uncontrolled areas, where controlled and uncontrolled areas are as defined in NCRP 147. Unless otherwise stated, the prescribed Pb for a wall is to span from the floor to 7 ft above the floor. Additionally, the Pb may be installed on either side of the wall.

No responsibility is assumed or implied for the adequacy of the enclosed specifications for workloads or conditions other than those stated.

Thank you for using Radiological Physics Associates for your physics consulting needs. If there are any questions concerning the design or if we can be of any further assistance, please feel free to contact me.

Respectfully,

A handwritten signature in black ink that reads "Alexander Woods".

Alexander Woods, MS, PhD  
Physicist



**Radiological Physics Associates**  
Medical Physics Consulting

2603 Elm Hill Pike, Suite P  
Nashville, TN 37214

800-262-6723  
inspections@rpanashville.com  
www.rpanashville.com

## Radiation Shielding Plan

**Pikeville Urgent Care**  
**3062 Main Street**  
**Pikeville, Tennessee 37367**

**January 2, 2026**

**DATE OF INSPECTION**

**Alexander Woods, MS, PhD**

**PHYSICIST**

**FACILITY**

**Toshiba Aquilion Prime**

**CT**

**EQUIPMENT**

**ROOM**

**Beam Orientation:** Tube directed toward gantry

**Radiographic Workload:** 180 Patients/Week

Barrier Description	Occupancy Type (C or UC)	Barrier Type (P or S)	Occupation Factor (T)	Use Factor (U)	Distance from Source (ft)	Maximum Allowable Exposure (mR/week)	Prescribed Shielding (inches Pb)
W1 - Control	C	S	1.00	1.00	10.0	<10	1/16" Pb
W2 - Exam	UC	S	0.50	1.00	7.0	<2	1/16" Pb
W3 - Corridor	UC	S	0.20	1.00	5.0	<2	1/16" Pb
D1 & D2 - Corridor	UC	S	0.20	1.00	5.0	<2	1/16" Pb
W4 - X-Ray	UC	S	1.00	1.00	6.0	<2	3/32" Pb
W5 - Exterior	UC	S	0.03	1.00	6.0	<2	*1/16" Pb*
Floor	UC	S	0.03	1.00	14.0	<2	N/A
Ceiling	UC	S	0.03	1.00	14.0	<2	N/A

**Ceilings and Floors:** No occupancy is assumed above or below the scan room.

**Comments:** All viewing windows must have same or greater mm Pb equivalence as corresponding barrier.  
\* Assuming brick or CMU is used for the exterior walls, only 1/16" lead is recommended at the windows. For any sections with less than 2.5" concrete or brick, 1/16" lead is recommended.



**Radiological Physics Associates**  
Medical Physics Consulting

2603 Elm Hill Pike, Suite P  
Nashville, TN 37214

800-262-6723   
inspections@rpanashville.com   
www.rpanashville.com 

---

Assumptions for Design and Construction

Windows within control booths for observing patients must provide the same degree of shielding as the barrier. The When adding sheet Pb, the Pb exhibits 'cold flow' if improperly secured. The Pb should be laminated to sheet rock or contact cemented to the wall. Do not simply nail the Pb to the walls. Overlap of the adjoining sheets be a half inch or more.

If poured concrete barriers are constructed, voids must be prevented. If using solid concrete bricks, take care that a solid course of mortar joins the individual bricks. Do not 'spot' cement the bricks.

Vents, pipes, or other violation of the shielding integrity should enter the room above the 7 ft height limitation. Standard wood doors may be used, if the Pb is laminated or contact cemented to the door. The Pb should be equivalent to the calculated barrier.

Corridors, toilets, dressing rooms, and other areas not directly related to the use of radiation are considered uncontrolled areas.




# Radiological Physics Associates

Medical Physics Consulting

2603 Elm Hill Pike, Suite P  
Nashville, TN 37214

800-262-6723  
inspections@rpanashville.com  
www.rpanashville.com





## Summary of Calculations

For the purposes of shielding calculations for CT scanners, only secondary radiation is considered, as the primary beam is well attenuated by the gantry. NCRP 147 demonstrates that the scattered air kerma per patient at 1 meter can be approximated as

$$K_S^1 = \kappa \frac{L}{p} CTDI_{100}$$

$L$  is the length scanned.  $p$  is known as the pitch and is defined as  $b/T_b$ .  $b$  is the patient movement through the bore per gantry rotation, and  $T_b$  is the beam width. For head and body phantoms, the scatter fraction per centimeter for the peripheral axis is assumed to be

$$\kappa_{\text{head}} = 9 \times 10^{-5} \text{ cm}^{-1}$$

$$\kappa_{\text{body}} = 3 \times 10^{-4} \text{ cm}^{-1}$$

For CT installations, it is convenient to use the dose length product,  $DLP$ , to compute shielding requirements, where

$$DLP = CTDI_{\text{vol}} L$$

and

$$CTDI_{\text{vol}} = \frac{1/3 CTDI_{100, \text{center}} + 2/3 CTDI_{100, \text{periphery}}}{p}$$

By assuming that the  $CTDI_{100}$  for the peripheral axis in the body phantom is twice the  $CTDI_{100}$  for the central axis and the  $CTDI_{100}$  for the peripheral axis and central axis are equal in the head phantom,

$$K_{\text{sec}}^1 (\text{head}) = \kappa_{\text{head}} DLP$$

$$K_{\text{sec}}^1 (\text{body}) = 1.2 \kappa_{\text{body}} DLP$$

Given a  $DLP$ , the air kerma can be scaled to appropriate distances and workloads to be included in a broad-beam transmission function. Thus, Figures A.2 and A.3 in Appendix A of NCRP 147 can be used to provide a barrier thickness to meet shielding design goals.



# Radiological Physics Associates

Medical Physics Consulting

2603 Elm Hill Pike, Suite P  
Nashville, TN 37214

800-262-6723  
inspections@rpanashville.com  
www.rpanashville.com

## Summary of Calculations

Barrier Description: W1 - Control

Source Distance: 3.0 m

Source Modality: CT

Occupancy Type (P):

Occupation Factor (T):

0.10 mGy per week

1.00

$$K_{sec}^{1}(\text{body}) = 0.28 \text{ mGy procedure}^{-1}$$

$$K_{sec}^{1}(\text{head}) = 0.15 \text{ mGy procedure}^{-1}$$

Material: Lead

$\alpha$ (mm <sup>-1</sup> )	2.246
$\beta$ (mm <sup>-1</sup> )	5.73
$\gamma$	0.547
K (mGy week <sup>-1</sup> )	5.01
B	0.01998
$x_{\text{barrier}}$ (mm)	0.924

$x_{\text{total}}$	0.924 mm	Thickness
$x_{\text{total}}$	2/55 inch	Required
$x_{\text{total}}$	1/16 inch	Recommended

Comments:



# Radiological Physics Associates

Medical Physics Consulting

2603 Elm Hill Pike, Suite P  
Nashville, TN 37214

800-262-6723  
inspections@rpanashville.com  
www.rpanashville.com

## Summary of Calculations

Barrier Description: W2 - Exam

Source Distance: 2.1 m

Source Modality: CT

Occupancy Type (P):

Occupation Factor (T):

0.02 mGy per week

0.50

$$K_{sec}^{1}(\text{body}) = 0.28 \text{ mGy procedure}^{-1}$$

$$K_{sec}^{1}(\text{head}) = 0.15 \text{ mGy procedure}^{-1}$$

Material: Lead

$\alpha$ (mm <sup>-1</sup> )	2.246
$\beta$ (mm <sup>-1</sup> )	5.73
$\gamma$	0.547
K (mGy week <sup>-1</sup> )	10.21
B	0.00392
$x_{\text{barrier}}$ (mm)	1.531

$x_{\text{total}}$	1.531 mm	Thickness
$x_{\text{total}}$	5/83 inch	Required
$x_{\text{total}}$	1/16 inch	Recommended

Comments:



# Radiological Physics Associates

Medical Physics Consulting

2603 Elm Hill Pike, Suite P  
Nashville, TN 37214

800-262-6723  
inspections@rpanashville.com  
www.rpanashville.com

## Summary of Calculations

Barrier Description: W3 - Corridor

Source Distance: 1.5 m

Source Modality: CT

Occupancy Type (P):

Occupation Factor (T):

0.02 mGy per week

0.20

$$K^1_{sec}(body) = 0.28 \text{ mGy procedure}^{-1}$$

$$K^1_{sec}(head) = 0.15 \text{ mGy procedure}^{-1}$$

Material: Lead

$\alpha$ (mm <sup>-1</sup> )	2.246
$\beta$ (mm <sup>-1</sup> )	5.73
$\gamma$	0.547
K (mGy week <sup>-1</sup> )	20.02
B	0.00499
$x_{barrier}$ (mm)	1.435

$x_{total}$	1.435 mm	Thickness
$x_{total}$	3/53 inch	Required
$x_{total}$	1/16 inch	Recommended

Comments:





# Radiological Physics Associates

Medical Physics Consulting

2603 Elm Hill Pike, Suite P  
Nashville, TN 37214

800-262-6723  
inspections@rpanashville.com  
www.rpanashville.com

## Summary of Calculations

Barrier Description: D1 & D2 - Corridor

Source Distance: 1.5 m  
Source Modality: CT

Occupancy Type (P): 0.02 mGy per week  
Occupation Factor (T): 0.20

$$K_{sec}^{1}(\text{body}) = 0.28 \text{ mGy procedure}^{-1}$$
$$K_{sec}^{1}(\text{head}) = 0.15 \text{ mGy procedure}^{-1}$$

Material: Lead

$\alpha$ (mm <sup>-1</sup> )	2.246
$\beta$ (mm <sup>-1</sup> )	5.73
$\gamma$	0.547
K (mGy week <sup>-1</sup> )	20.02
B	0.00499
$x_{\text{barrier}}$ (mm)	1.435

$x_{\text{total}}$	1.435 mm	Thickness
$x_{\text{total}}$	3/53 inch	Required
$x_{\text{total}}$	1/16 inch	Recommended

Comments:



# Radiological Physics Associates

Medical Physics Consulting

2603 Elm Hill Pike, Suite P  
Nashville, TN 37214

800-262-6723  
inspections@rpanashville.com  
www.rpanashville.com

## Summary of Calculations

Barrier Description: W4 - X-Ray

Source Distance: 1.8 m

Source Modality: CT

Occupancy Type (P):

Occupation Factor (T):

0.02 mGy per week

1.00

$$K_{sec}^{1} (body) = 0.28 \text{ mGy procedure}^{-1}$$

$$K_{sec}^{1} (head) = 0.15 \text{ mGy procedure}^{-1}$$

Material: Lead

$\alpha$ (mm <sup>-1</sup> )	2.246
$\beta$ (mm <sup>-1</sup> )	5.73
$\gamma$	0.547
K (mGy week <sup>-1</sup> )	13.90
B	0.00144
$x_{barrier}$ (mm)	1.938

$x_{total}$	1.938 mm	Thickness
$x_{total}$	1/13 inch	Required
$x_{total}$	3/32 inch	Recommended

Comments:

1/16" lead-lined wallboard at each side of the wall (CT & X-Ray) is recommended.



# Radiological Physics Associates

Medical Physics Consulting

2603 Elm Hill Pike, Suite P  
Nashville, TN 37214

800-262-6723  
inspections@rpanashville.com  
www.rpanashville.com

## Summary of Calculations

Barrier Description: W5 - Exterior

Source Distance: **1.8 m**

Source Modality: **CT**

Occupancy Type (P):

**0.02 mGy per week**

Occupation Factor (T):

**0.03**

$$K_{sec}^{1}(body) = 0.28 \text{ mGy procedure}^{-1}$$

$$K_{sec}^{1}(head) = 0.15 \text{ mGy procedure}^{-1}$$

Material: Concrete

$\alpha$ (mm <sup>-1</sup> )	0.0383
$\beta$ (mm <sup>-1</sup> )	0.0142
$\gamma$	0.658
K (mGy week <sup>-1</sup> )	13.90
B	0.05754
<b><math>x_{barrier}</math> (mm)</b>	<b>64.223</b>

<b><math>x_{total}</math></b>	<b>64.223 mm</b>	Thickness
<b><math>x_{total}</math></b>	<b>2 28/53 inch</b>	Required
<b><math>x_{total}</math></b>	<b>2.5 inch</b>	<b>Recommended</b>

Comments:



# Radiological Physics Associates

Medical Physics Consulting

2603 Elm Hill Pike, Suite P  
Nashville, TN 37214

800-262-6723  
inspections@rpanashville.com  
www.rpanashville.com

## Summary of Calculations

Barrier Description: Floor

Source Distance: 4.3 m

Source Modality: CT

Occupancy Type (P):

0.02 mGy per week

Occupation Factor (T):

0.03

$$K_{sec}^{1}(\text{body}) = 0.28 \text{ mGy procedure}^{-1}$$

$$K_{sec}^{1}(\text{head}) = 0.15 \text{ mGy procedure}^{-1}$$

Material: Concrete

$\alpha$ (mm <sup>-1</sup> )	0.0383
$\beta$ (mm <sup>-1</sup> )	0.0142
$\gamma$	0.658
K (mGy week <sup>-1</sup> )	2.55
B	0.31327
$x_{\text{barrier}}$ (mm)	24.114

$x_{\text{total}}$	24.114 mm	Thickness
$x_{\text{total}}$	75/79 inch	Required
$x_{\text{total}}$	0.9 inch	Recommended

Current Thickness of Concrete in Floor:

Additional Required Thickness of Floor:

Comments: *There is no occupancy below. No additional shielding is recommended.*



# Radiological Physics Associates

Medical Physics Consulting

2603 Elm Hill Pike, Suite P  
Nashville, TN 37214

800-262-6723  
inspections@rpanashville.com  
www.rpanashville.com

## Summary of Calculations

Barrier Description: Ceiling

Source Distance: 4.3 m

Source Modality: CT

Occupancy Type (P):

0.02 mGy per week

Occupation Factor (T):

0.03

$$K_{sec}^{1}(body) = 0.28 \text{ mGy procedure}^{-1}$$

$$K_{sec}^{1}(head) = 0.15 \text{ mGy procedure}^{-1}$$

Material: Concrete

$\alpha$ (mm <sup>-1</sup> )	0.0383
$\beta$ (mm <sup>-1</sup> )	0.0142
$\gamma$	0.658
K (mGy week <sup>-1</sup> )	2.55
B	0.31327
$x_{barrier}$ (mm)	24.114

$x_{total}$	24.114 mm	Thickness
$x_{total}$	75/79 inch	Required
$x_{total}$	0.9 inch	Recommended

Current Thickness of Concrete in Ceiling:

Additional Required Thickness of Ceiling:

Comments: *There is no occupancy above. No additional shielding is recommended.*